

goss、cube、copper 30%を元に

各種 α 範囲の極点図からTexToolsとLaboTexのVolumeFractionを計算した。
 LaboTexのVとFは広がりを変可変と10degに固定を示す。
 固定はTexTools,LaboTex共に、10degとした。

	作成	goss	cube	copper
再計算極点図	LaboTexVF-V	28.70	29.60	31.50
	LaboTexVF-F	24.87	26.54	24.64
	TexTools	23.46	23.93	23.66
0->75	LaboTexVF-V	28.50	29.40	32.10
	LaboTexVF-F	24.79	26.52	24.60
	TexTools	16.05	16.06	16.00
15->90	LaboTexVF-V	30.00	29.50	30.50
	LaboTexVF-F	24.86	26.56	24.66
	TexTools	23.50	24.15	23.67
15->75	LaboTexVF-V	28.60	29.30	31.90
	LaboTexVF-F	24.76	26.51	24.63
	TexTools	16.06	16.28	16.03
{111}0->90	LaboTexVF-V	29.10	28.90	31.90
	LaboTexVF-F	24.78	26.56	24.55
	TexTools	23.37	24.14	23.42
{111}0->75	LaboTexVF-V	29.00	29.10	30.80
	LaboTexVF-F	23.45	26.67	24.36
	TexTools	16.03	15.95	15.37
{111}15->75	LaboTexVF-V	29.00	29.10	30.80
	LaboTexVF-F	23.40	26.53	24.29
	TexTools	15.93	15.84	15.19


TexToolsは α 範囲を制限すると、VolumeFractionが減少する。
 VolumeFractionが減少しているのではなく、広がりが大きくなっている

この現象で、実際の測定結果から、TexToolsはODFが滑らかになる事が理解できる。


2010年03月04日

HelperTex


goss, cube, copper 方位を各々30%としてODF計算



No.	VF(%)	Phi1(FWHM)	Phi(FWHM)	Phi2(FWHM)	Orientation
1:	30.0	10.0	10.0	10.0	{ 1 1 0 } _K 0 0 1 > goss
2:	30.0	10.0	10.0	10.0	{ 0 0 1 } _K 1 0 0 > cube
3:	30.0	10.0	10.0	10.0	{ 1 1 2 } _K 1 1 -1 > copper
4:	10.01	Background Volume Fraction			



{111}、{220}、{311} 再計算極点図を作成




再計算極点図を export し再度ODFに読み込んでVolume Fractionを計算

SelectFile(TXT(b,intens),TXT2(a,b,intens.))	h,k,l	2Theta	Alfa Area	AlfaS	AlfaE	Select
111_5.0deg-rp2.TXT	1,1,1	0.0	0.0->90.0	0.0	90.0	<input checked="" type="checkbox"/>
220_5.0deg-rp2.TXT	2,2,0	0.0	0.0->90.0	0.0	90.0	<input checked="" type="checkbox"/>
311_5.0deg-rp2.TXT	3,1,1	0.0	0.0->90.0	0.0	90.0	<input checked="" type="checkbox"/>

α 範囲を 0 → 90 の場合


LaboTexでは



No.	VF(%)	Phi1(FWHM)	Phi(FWHM)	Phi2(FWHM)	Orientation
1:	28.7	10.4	9.6	9.6	{ 1 1 0 } < 0 0 1 > goss
2:	29.6	10.0	9.9	10.3	{ 0 0 1 } < 1 0 0 > cube
3:	31.5	11.7	9.3	10.1	{ 1 1 2 } < 1 1 -1 > copper
4:	10.20	Background Volume Fraction			

Volume Fraction	Delta Phi1	Delta Phi	Delta Phi2	Orientation
24.87	10.00	10.00	10.00	{ 1 1 0 } < 0 0 1 > goss
26.54	10.00	10.00	10.00	{ 0 0 1 } < 1 0 0 > cube
24.64	10.00	10.00	10.00	{ 1 1 2 } < 1 1 -1 > copper

TextToolsでは



Orientation Volume Fraction

Ideal Orientations

All of above

ϕ_1 # Φ # ϕ_2 #

Tolerance (degree) 10

With sample symmetrical orientations (4 fold)

Fibers

gamma fiber (phi = 55, phi2 = 45)

h 1 k 1 l 1

Tolerance (degree) 15

ODF File Name G:\DATA\GCC30\3peak.HODF >>


Result: Cube=23.46%Goss=23.93%Brass=0.09%S=1.27%Copper=23.66%R-Cube=0.05%

Calculate

Cancel

α 範囲 0 → 7.5 度の場合


SelectFile(TXT(b,intens),TXT2(a,b,intens))	h,k,l	2Theta	Alfa Area	AlfaS	AlfaE	Select
 111_5.0deg-rp2.TXT	1,1,1	0.0	0.0→90.0	0.0	75	<input checked="" type="checkbox"/>
 220_5.0deg-rp2.TXT	2,2,0	0.0	0.0→90.0	0.0	75	<input checked="" type="checkbox"/>
 311_5.0deg-rp2.TXT	3,1,1	0.0	0.0→90.0	0.0	75	<input checked="" type="checkbox"/>



No.	VF(%)	Phi1(FWHM)	Phi(FWHM)	Phi2(FWHM)	Orientation
1:	28.5	10.5	9.5	9.5	{ 1 1 0 } < 0 0 1 > goss
2:	29.4	10.0	9.8	10.2	{ 0 0 1 } < 1 0 0 > cube
3:	32.1	11.9	9.1	10.2	{ 1 1 2 } < 1 1 -1 > copper
4:	10.06	Background Volume Fraction			

Volume Fraction	Delta Phi1	Delta Phi	Delta Phi2	Orientation
24.79	10.00	10.00	10.00	{ 1 1 0 } < 0 0 1 > goss
26.52	10.00	10.00	10.00	{ 0 0 1 } < 1 0 0 > cube
24.60	10.00	10.00	10.00	{ 1 1 2 } < 1 1 -1 > copper

TextToolsでは



Orientation Volume Fraction

Ideal Orientations

All of above

ϕ_1 # Φ # ϕ_2 #

Tolerance (degree) 10

With sample symmetrical orientations (4 fold)

Fibers

gamma fiber (phi = 55, phi2 = 45)

h | k | l | l

Tolerance (degree) 15

ODF File Name G:\DATA\GCC30\0-75.HODF >>




Result: Cube=16.05%Goss=16.06%Brass=0.10%S=2.31%Copper=16.00%R-Cube=0.05%

Calculate

Cancel

α 範囲を 15 → 90 の場合


SelectFile(TXT(b,intens),TXT2(a,b,intens.))	h,k,l	2Theta	Alfa Area	AlfaS	AlfaE	Select
 111_5.0deg-rp2.TXT	1,1,1	0.0	0.0->90.0	15	90.0	<input checked="" type="checkbox"/>
 220_5.0deg-rp2.TXT	2,2,0	0.0	0.0->90.0	15	90.0	<input checked="" type="checkbox"/>
 311_5.0deg-rp2.TXT	3,1,1	0.0	0.0->90.0	15	90.0	<input checked="" type="checkbox"/>



No.	VF(%)	Phi1(FWHM)	Phi(FWHM)	Phi2(FWHM)	Orientation
1:	30.0	10.9	9.6	9.5	{ 1 1 0 } < 0 0 1 > goss
2:	29.5	10.0	9.8	10.3	{ 0 0 1 } < 1 0 0 > cube
3:	30.5	11.1	9.3	9.9	{ 1 1 2 } < 1 1 -1 > copper
4:	10.08	Background Volume Fraction			

Volume Fraction	Delta Phi1	Delta Phi	Delta Phi2	Orientation
24.86	10.00	10.00	10.00	{ 1 1 0 } < 0 0 1 > goss
26.56	10.00	10.00	10.00	{ 0 0 1 } < 1 0 0 > cube
24.66	10.00	10.00	10.00	{ 1 1 2 } < 1 1 -1 > copper

TextToolsでは



Orientation Volume Fraction

Ideal Orientations

All of above

ϕ_1 # Φ # ϕ_2 #

Tolerance (degree) 10

With sample symmetrical orientations (4 fold)

Fibers

gamma fiber (phi = 55, phi2 = 45)

h | k | l |

Tolerance (degree) 15

ODF File Name G:\DATA\GCC30%\15-90\15-90.HODF >>


Result: Cube=23.50%Goss=24.15%Brass=0.09%S=1.28%Copper=23.67%R-Cube=0.05%

Calculate

Cancel

α 範囲 1 5 - 7 5 の場合


SelectFile(TXT(b,intens),TXT2(a,b,intens))	h,k,l	2Theta	Alfa Area	AlfaS	AlfaE	Select
 111_5.0deg-rp2.TXT	1,1,1	0.0	0.0->90.0	15	75	<input checked="" type="checkbox"/>
 220_5.0deg-rp2.TXT	2,2,0	0.0	0.0->90.0	15	75	<input checked="" type="checkbox"/>
 311_5.0deg-rp2.TXT	3,1,1	0.0	0.0->90.0	15	75	<input checked="" type="checkbox"/>



No.	VF(%)	Phi1(FWHM)	Phi(FWHM)	Phi2(FWHM)	Orientation
1:	28.6	10.7	9.6	9.5	{ 1 1 0 } < 0 0 1 > goss
2:	29.3	9.9	9.8	10.2	{ 0 0 1 } < 1 0 0 > cube
3:	31.9	12.0	9.2	10.1	{ 1 1 2 } < 1 1 -1 > copper
4:	10.17	Background Volume Fraction			

Volume Fraction	Delta Phi1	Delta Phi	Delta Phi2	Orientation
24.76	10.00	10.00	10.00	{ 1 1 0 } < 0 0 1 > goss
26.51	10.00	10.00	10.00	{ 0 0 1 } < 1 0 0 > cube
24.63	10.00	10.00	10.00	{ 1 1 2 } < 1 1 -1 > copper

TextTools



Orientation Volume Fraction

Ideal Orientations

All of above

ϕ_1 # Φ # ϕ_2 #

Tolerance (degree) 10

With sample symmetrical orientations (4 fold)

Fibers

gamma fiber (phi = 55, phi2 = 45)

h | k | l | l

Tolerance (degree) 15

ODF File Name G:\DATA\GCC30\15-90\15-75.HODF >>




Result: Cube=16.06%Goss=16.28%Brass=0.08%S=2.32%Copper=16.03%R-Cube=0.05%

Calculate

Cancel

{ 1 1 1 } のみの場合


SelectFile(TXT(b,intens),TXT2(a,b,intens.))	h,k,l	2Theta	Alfa Area	AlfaS	AlfaE	Select
 111_50deg-rp2.TXT	1,1,1	0.0	0.0->90.0	0.0	90.0	<input checked="" type="checkbox"/>
	2,0,0	0.0		0.0	60.0	<input type="checkbox"/>



No.	VF(%)	Phi1(FWHM)	Phi(FWHM)	Phi2(FWHM)	Orientation
1:	29.1	11.3	9.4	9.5	{ 1 1 0 } < 0 0 1 > goss
2:	28.9	9.5	10.3	10.2	{ 0 0 1 } < 1 0 0 > cube
3:	31.9	11.5	9.4	10.4	{ 1 1 2 } < 1 1 -1 > copper
4:	10.12	Background Volume Fraction			

Volume Fraction	Delta Phi1	Delta Phi	Delta Phi2	Orientation
24.78	10.00	10.00	10.00	{ 1 1 0 } < 0 0 1 > goss
26.56	10.00	10.00	10.00	{ 0 0 1 } < 1 0 0 > cube
24.55	10.00	10.00	10.00	{ 1 1 2 } < 1 1 -1 > copper

TextTools



Orientation Volume Fraction ✖

Ideal Orientations Fibers

All of above gamma fiber (phi = 55, phi2 = 45)

ϕ_1 # Φ # ϕ_2 #

h k l

Tolerance (degree) 10 Tolerance (degree) 15

With sample symmetrical orientations (4 fold)


ODF File Name G:\DATA\GCC30\0-90-111.HODF >>

Result: Cube=23.37%Goss=24.14%Brass=0.11%S=1.39%Copper=23.42%R-Cube=0.04%

Calculate
Cancel

α 範圍 0 - 7.5


SelectFile(TXT(b,intens),TXT2(a,b,intens.))	h,k,l	2Theta	Alfa Area	AlfaS	AlfaE	Select
 111_5.0deg-rp2.TXT	1,1,1	0.0	0.0->90.0	0.0	75	<input checked="" type="checkbox"/>
	2,0,0	0.0		0.0	60.0	<input type="checkbox"/>



No.	VF(%)	Phi1(FWHM)	Phi(FWHM)	Phi2(FWHM)	Orientation
1:	29.0	14.0	9.0	9.5	{ 1 1 0 } < 0 0 1 > goss
2:	29.1	9.7	10.4	10.3	{ 0 0 1 } < 1 0 0 > cube
3:	30.8	11.4	9.5	9.9	{ 1 1 2 } < 1 1 -1 > copper
4:	0.0	12.6	12.4	13.7	{ 0 0 1 } < 5 1 0 >
5:	11.12	Background Volume Fraction			

Volume Fraction	Delta Phi1	Delta Phi	Delta Phi2	Orientation
23.45	10.00	10.00	10.00	{ 1 1 0 } < 0 0 1 > goss
26.67	10.00	10.00	10.00	{ 0 0 1 } < 1 0 0 > cube
24.36	10.00	10.00	10.00	{ 1 1 2 } < 1 1 -1 > copper

TextTools



Orientation Volume Fraction

Ideal Orientations Fibers

All of above

ϕ_1 # Φ # ϕ_2 #

Tolerance (degree) 10

With sample symmetrical orientations (4 fold)

gama fiber (phi = 55, phi2 = 45)

h | k | l |

Tolerance (degree) 15

ODF File Name G:\DATA\GCC30\0-75-111\0-75-111.HODF >>

Result: Cube=16.03%Goss=15.95%Brass=0.08%S=2.23%Copper=15.37%R-Cube=0.04%


Calculate

Cancel

The software interface includes a 'Result' field displaying the following data: Cube=16.03%, Goss=15.95%, Brass=0.08%, S=2.23%, Copper=15.37%, R-Cube=0.04%. Below the result field is a horizontal bar chart consisting of 15 blue rectangular segments of varying heights, representing the distribution of the different orientation types.

α 範圍 1 5 - 7 5 - { 1 1 1 }


SelectFile(TXT(b,intens),TXT2(a,b,intens.))	h,k,l	2Theta	Alfa Area	AlfaS	AlfaE	Select
 111_5.0deg-rp2.TXT	1,1,1	0.0	0.0->90.0	15	75	<input checked="" type="checkbox"/>
	2,0,0	0.0		0.0	60.0	<input type="checkbox"/>



No.	VF(%)	Phi1(FWHM)	Phi(FWHM)	Phi2(FWHM)	Orientation
1:	29.0	14.0	9.0	9.5	{ 1 1 0 } < 0 0 1 > goss
2:	29.1	9.7	10.4	10.3	{ 0 0 1 } < 1 0 0 > cube
3:	30.8	11.4	9.5	9.9	{ 1 1 2 } < 1 1 -1 > copper
4:	0.0	12.6	12.4	13.7	{ 0 0 1 } < 5 1 0 >
5:	11.12	Background Volume Fraction			

Volume Fraction	Delta Phi1	Delta Phi	Delta Phi2	Orientation
23.40	10.00	10.00	10.00	{ 1 1 0 } < 0 0 1 > goss
26.53	10.00	10.00	10.00	{ 0 0 1 } < 1 0 0 > cube
24.29	10.00	10.00	10.00	{ 1 1 2 } < 1 1 -1 > copper

TextTools



Orientation Volume Fraction ✖

Ideal Orientations Fibers

All of above gamma fiber (phi = 55, phi2 = 45)

ϕ_1 # Φ # ϕ_2 #

h k l

Tolerance (degree) 10 Tolerance (degree) 15

With sample symmetrical orientations (4 fold)

ODF File Name G:\DATA\GCC30\15-75\15-75.HODF >>

Result: Cube=15.93%Goss=15.84%Brass=0.08%S=2.25%Copper=15.19%R-Cube=0.04%

Calculate
Cancel