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


The screenshot shows the ODFPolefigure1.5 software interface. The 'Background delete mode' section has 'DoubleMod...' checked. The 'Smoothing' section has 'Full' checked. The 'AbsCalc' section has 'Schulz reflection method' selected. The 'Defocus file' section has 'Normalization' checked. The 'Smoothing(for ADC)' section has 'Cycles' set to 2 and 'Weight' set to 15. The 'Normalization' section has 'CTR' checked. The 'CenterData' section has 'Average' checked. The 'OutFiles' section has 'TXT2' selected. The 'ValueODF-B' and 'ValueODF-A' buttons are visible. The bottom right corner shows 'Select crystal : Cubic' and '19/11/14'.

## ODFPoleFigure2 処理

ODFPoleFigure2 3.92 by CTR PDuser ctr helpertex

File Linear(absolute)3D ToolKit Help InitSet BGMode Measure Condition Free OverlapRevision MinimumMode Rp% Normalization

Files select: ASC(RINT-PC) 110\_TEST\_θ group-chback4.asc

Calculation Condition: C:\CTR\tmp\group\110\_TEST\_θ group-chback4.asc

Background delete mode:  DoubleMode  SingleMode  LowMode  HighMode  Nothing

BG defocus: DSH1.2mm+Schulz+RSH5mm

Smoothing:  +α  3  Arithmetic mean

Peak slit: 2.0 mm BG Slit: 2.0 mm  Peak Slit / BGS.. BG Scope: 80.0 deg. 90.0 deg.

AbsCalc:  Ref  Trans Schulz reflection method Absorption coefficient: 13.9 1/cm Thickness: 0.1 cm

Defocus file Select Transmission defocus HKL+T


Defocus(3) function files folder(Calc unbackdefocus): BB185mm Limit Alfa Defocus value: Free(LimitValue=0.0)

Smoothing for ADC: Cycles 2 Weight 10

Normalization:  CTR  Connect  Average

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