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ODFPoleFigure2 3.95T[20/10/31] by CTR

File Linear(absolute)Contour ToolKit Help InitSet BGMode Measure Condition Free OverlapRevision MinimumMode Rp% Normalization

Files select: ASC(RINT-PC) 110.ASC 200.ASC 211.ASC 222.ASC

Calculation Condition: C:\CTR\DATA\TEST\PCDATA\random#110.ASC

Background delete mode:  DoubleMode  SingleMode  LowMode  HighMode  Nothing

BG defocus: DSH1.2mm+Schulz+RSH5mm

Smoothing:  +α 3 Arithmetic mean Disp

RD:  0.0 Interpolation Full Disp

Peak slit: 7.0 mm BG Slit: 7.0 mm  PeakSlit / BGS\_ BG Scope: 80.0 deg. 90.0 deg. Set Disp  Inhibit

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Defocus file Select Transmission defocus HKL+T

Defocus(1) functions file C:\CTR\DATA\TEST\PCDATA\random\defocus\DEFOCUS\_NOTNORM\_F.TXT

Make defocus function files by TXT2 Files  Normalization degree of a polynomial 0  TenckhoffFitting

Defocus(3) function files folder(Calc unbackdefocus) BB185mm Limit Alfa Defocus value Free(LimitValue=0.0)

Defocus(2) function files folder(Calc backdefocus) DSH1.2mm+Schulz+RSH5mm  Search minimum EqualAngleRp%(Cubic only)

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The screenshot displays the ODF software interface with four circular plots at the top showing intensity distributions for different hkl values: (1,1,0), (2,0,0), (2,1,1), and (2,2,2). The main window is titled "ODFPoleFigure2 3.95T[20/10/31] by CTR". The menu bar includes File, Linear(absolute)Contour, Toolkit, Help, InitSet, BGMode, Measure, Condition, Free, OverlapRevision, MinimumMode, Rp%, and Normalization. The "Files select" section shows "ASC(RINT-PC)" and a list of files: 110.ASC, 200.ASC, 211.ASC, 222.ASC. The "Calculation Condition" section shows the path "C:\CTR\DATA\TEST\PCDATA\sample#110.ASC" and "hkl 1,1,0". The "Background delete mode" section has "LowMode" selected. The "Smoothing" section has "Arithmetic mean" selected. The "AbsCalc" section has "Schulz reflection method" selected. The "Defocus file" section has "Defocus(1) functions file" selected. The "Normalization" section has "CTR" selected. The status bar shows "Filemake success !!", "CTRHome : C:", "Select crystal : Cubic", and "20/08/26".

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The screenshot displays the ODF software interface with four circular plots at the top showing intensity distributions for different hkl values: (1,1,0), (2,0,0), (2,1,1), and (2,2,2). The main window is titled "ODFPoleFigure2 3.95T[20/10/31] by CTR". The menu bar includes File, Linear(absolute)Contour, Toolkit, Help, InitSet, BGMode, Measure, Condition, Free, OverlapRevision, MinimumMode, Rp%, and Normalization. The "Files select" section shows "ASC(RINT-PC)" and a list of files: 110.ASC, 200.ASC, 211.ASC, 222.ASC. The "Calculation Condition" section shows the path "C:\CTR\DATA\TEST\PCDATA\sample#110.ASC" and "hkl 1,1,0". The "Background delete mode" section has "LowMode" selected. The "Smoothing" section has "Arithmetic mean" selected. The "AbsCalc" section has "Schulz reflection method" selected. The "Defocus file" section has "Defocus(1) functions file" selected. The "Normalization" section has "CTR" selected. The status bar shows "Filemake success !!", "CTRHome : C:", "Select crystal : Cubic", and "20/08/26".

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ODFPoleFigure2 3.95T[20/10/31] by CTR

File Linear(absolute)Contour Toolkit Help InitSet BGMode Measure Condition Free OverlapRevision MinimumMode Rp% Normalization

Files select: ASC(RINT-PC) 110.ASC 200.ASC 211.ASC 222.ASC

Calculation Condition: Previous Next C:\CTR\DATA\TEST\PCDATA\sample#110.ASC hkl 1,1,0

Background delete mode:  DoubleMode  SingleMode  LowMode  HighMode  Nothing BG defocus DSH1.2mm+Schulz+RSH5mm

Minimum( $\alpha, \beta$ )  MinimumAverage( $\alpha, \beta$ )  0.5 Trans blinds angle 30.0

Peak slit 7.0 mm BG Slit 7.0 mm  PeakSlit / BGS.. BG Scope 80.0 deg. 90.0 deg.

AbsCalc:  Ref  Trans Schulz reflection method Absorption coefficient 303.149 1/cm Thickness 0.02 cm 2Theta 20.04 deg.  1/Kt Profile

Defocus file Select Transmission defocus HKL+T

Defocus(1) functions file C:\CTR\DATA\TEST\PCDATA\random#defocus\DEFOCUS\_NOTNORM\_F.TXT


Defocus(3) function files folder(Calc unbackdefocus) BB185mm Limit Alfa Defocus value Free(LimitValue=0.0)

Defocus(2) function files folder(Calc backdefocus) DSH1.2mm+Schulz+RSH5mm  Search minimum EqualAngleRp%(Cubic only)  1/Ra Profile



Smoothing for ADC:  Cycles 2 Weight 10

Normalization:  CTR  Connect  Average OutFiles:  Asc  MTextAsc  Ras  TXT  TXT2

Filemake success !! CTRHome : C: Select crystal : Cubic 20/08/26

Before

110	100	211	111
110	100	211	111
96.8	16.2	18.0	28.9


Average= 39.9 %

110	100	211	111
110	100	211	111
9.6	4.7	5.6	5.9

Average= 6.4 %

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PfToODF3 8.47T[20/10/31] by CTR

File Option Symmetric Software Data Help

Lattice constant

Material: A-Iron-Measure-IntegralData.txt

Structure Code(Symmetries after Schoenfiles): 7 - O (cubic)

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ODFPoleFigure2.35T[20/10/31] by CTR

File Linear(absolute)Contour ToolKit Help InitSet BGMode Measure Condition Free OverlapRevision MinimumMode Rp% Normalization

Files select  
ASC(RINT-FC) 110.ASC 200.ASC 211.ASC 222.ASC

Calculation Condition  
C:\CTR\DATA\TEST\PCDATA\sample\notrandom\110.ASC

Background delete mode  
 DoubleMode  SingleMode  LowMode  HighMode  Nothing  
 Minimum( $\alpha, \beta$ )  MinimumAverage( $\alpha, \beta$ ) 0.5  
 Peak slit 7.0 mm BG Slit 7.0 mm  PeakSlit / BGS... BG Scope 80.0 deg. 90.0 deg. Set Disp  $\alpha$ Inhibit

Smoothing  
 + $\alpha$  3 Arithmetic mean Disp  
 RD 0.0 Interpolation Full Disp

AbsCalc  
 Ref  Trans Schulz reflection method Change Absorption coefficient 303.149 1/cm Thickness 0.02 cm Set 2Theta 20.04 deg.  1/Kt Profile

Defocus file Select Transmission defocus HKL+T  
 Defocus(1) functions file C:\CTR\DATA\TEST\PCDATA\random\defocus\DEFOCUS.NOTFORM.F.TXT  
 Make defocus function files by TXT2 Files  Normalization degree of a polynomial 0  TenckhoffFitting TextDisp

Defocus(3) function files folder(Calc unbackdefocus) BB185mm Limit Alpha Defocus value Free(LimitValue=0.0)




Defocus(2) function files folder(Calc backdefocus) DSH1.2mm+Schulz+RSH5mm  Search minimum EqualAngleRp%(Cubic only)  1/Ra Profile

Smoothing for ADC  
 Cycles 2 Weight 10 Disp  
 After connection

Normalization  
 CTR  Connect  Average OutFiles  
 Asc  MTextAsc  Ras  TXT  TXT2

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