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

ODFPoleFigure2S 4.08T[23/12/31] by CTR

File Linear(absolute)3D Toolkit Help InitSet BGMode Defocus Condition Free OverlapRevision MinimumMode Rp% Normalization

Files select
ASC(RINT-PC) 111.ASC 200.ASC 220.ASC

Holder
Holder
Rename
Renamw-delete

Calculation Condition
Previous Next C:\CTR\DATA\Aluminum-H-O\Aluminum-O\111.ASC hkl 1,1,1 Change

Background delete mode
 DoubleMo... SingleMode LowMode HighMode Nothing BG defocus DSH1.2mm+Schulz+RSH5mm Minimum mo...
 Minimum($\alpha\beta$) MinimumAverage(α)X 0.5 Trans blinds angle 30.0
Peak slit 7.0 mm BG Slit 7.0 mm PeakSlit / BGS... BG Scope 80.0 deg. 90.0 deg. Set SetAll Disp α Inhibit

Smoothing
 + α 3 Arithmetic mean Disp
RD 0.0 Interpolation Full **Disp**

AbsCalc
 Ref Trans Schulz reflection method Change Absorption coefficient 133.0 1/cm Thickness 0.2 cm Set 2Theta 38.59 deg. 1/Kt Profile


Defocus file Select Transmission defocus HKL+T
 Defocus(1) functions file C:\tmp\Al-random\def\defocus\DEFOCUS_F.TXT
Make defocus function files by TXT2 Files Normalization degree of a polynomial 0 TenckhoffFitting TextDisp
 Defocus(3) function files folder(Calc unbackdefocus) BB285mm-Slit-1mm-7mm Limit Alfa Defocus value Free(LimitValue=0.0)
 Defocus(2) function files folder(Calc backdefocus) DSH1.2mm+Schulz+RSH5mm Search minimum EqualAngleRp%(Cubic only) 1/Ra Profile

Smoothing for ADC
 Cycles 2 Weight 10 Disp
 After connection

Normalization CenterData OutFiles
 CTR Connect Average ASC MTextAsc Ras TXT TXT2
Cancel Calc Connect
ValueODFVF-B ValueODFVF-A Average Exit&ODF ODF ALLODFFiles

CTRHome : C: Select crystal : Cubic 23/04/13

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The screenshot displays the ODFPoleFigure2S software interface. At the top, three circular pole figures are shown, followed by a 3D surface plot. The main window contains several panels and controls:

- File select:** Shows the current file path and holder information.
- Holder:** Includes buttons for Holder, Rename, and Renamw-delete.
- Calculation Condition:** Contains settings for Background delete mode (DoubleMo... selected), BG defocus (DSH1.2mm+Schulz+RSH5mm), and Smoothing (Arithmetic mean, Disp).
- AbsCalc:** Includes settings for Ref, Trans, Schulz reflection method, Absorption coefficient (133.0), Thickness (0.2), and 2Theta (38.59).
- Defocus file:** Includes settings for Defocus(1) functions file, Make defocus function files by TXT2, and Normalization (checked).
- Defocus(3) function files folder:** Includes settings for folder selection and Limit Alfa Defocus value.
- Defocus(2) function files folder:** Includes settings for folder selection and Search minimum EqualAngleRp% (checked).
- Smoothing for ADC:** Includes settings for Cycles (2), Weight (10), and Disp.
- Normalization:** Includes settings for CenterData (Average) and OutFiles (TXT2 selected).
- Buttons:** Includes Calc, Connect, and other control buttons.

Annotations in blue text point to specific settings:

- Background削除** points to the Background delete mode section.
- BG-defocus mde** points to the BG defocus dropdown menu.
- RD補正** points to the RD (0.0) field.
- Defocus補正** points to the Defocus(1) functions file section.
- 最適化Rp%** points to the Search minimum EqualAngleRp% checkbox.

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