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


The screenshot displays the ODF software interface with four circular diffraction patterns and a 3D surface plot. The interface includes various settings for background removal, defocus correction, and peak optimization. Key settings and annotations are as follows:

- Files select:** 111-NO 010.ASC 200-NO 010.ASC 220-NO 010.ASC 311-NO 010.ASC
- Holder:** Holder, Rename, Renamw-delete
- Calcration Condition:** Previous, Next, C:\TMP\AP\NO 010\ASC\220-NO 010.ASC
- Background delete mode:**  DoubleMode,  SingleMode,  LowMode,  HighMode,  Nothing
- BG defocus:** DSH1.2mm+Schulz+RSH5mm
- Smoothing:**  +a, 3, Arithmetic mean, Disp
- RD:**  0, Interpolation, Full, Disp
- AbsCalc:**  Ref,  Trans, Schulz reflection method, Absorption coefficient: 13.9, 1/cm, Thickness: 0.1, cm, 2Theta: 64.86, deg., 1/Kt, Profile
- Defocus file:**  Defocus(1) functions file, C:\TMP\AP\AI-powder-random\defocus\DEFocus\_F.TXT
- Make defocus function files by TXT2:** Files,  Normalization, degree of a polynomial: 0,  TenckhoffFitting, TextDisp
- Defocus(3) function files folder:** Calc unbackdefocus, BB185mm, Limit Alfa Defocus value: Free(LimitValue=0.0)
- Defocus(2) function files folder:** Calc backdefocus, DSH1.2mm+Schulz+RSH5mm,  Search minimum EqualAngleRp%(Cubic only), 1/Ra, Profile
- Smoothing for ADC:**  Cycles, 2, Weight: 10, Disp
- Normalization:**  CTR,  Connect,  Average
- CenterData:**  Average
- OutFiles:**  ASC,  MTextAsc,  Ras,  TXT,  TXT2
- Buttons:** Cancel, Calc, Connect
- Footer:** CTRHome : C:, Select crystal : Cubic, 23/04/14

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