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ODFPoleFigure2S 4.08 by CTR PDuser HelperTex CTR

File Linear(absolute)3D Toolkit Help InitSet BGMode Measure Condition Free OverlapRevision MinimumMode Rp% Normalization

Files select
ASC(RINT-PC) 111-7mm.ASC 200-7mm.ASC 220-7mm.ASC 311-7mm.ASC

Holder
Holder: C:\A1100P\N0004
Rename: C:\A1100P\N0005
Renamw-delete: C:\A1100P\N0006

Calcration Condition
Previous Next C:\A1100P\N0001\111-7mm.ASC

Background delete mode
 DoubleMo... SingleMode LowMode HighMode Nothing BG defocus: DSH1.2mm+Schulz+RSH5mm

Smoothing
 + α 3 Arithmetic mean Disp

RD
 0.0 Interporation Full Disp

Peak slit 7.0 mm BG Slit 7.0 mm PeakSlit / BGS... BG Scope 80.0 deg. 90.0 deg. Set SetAll Disp α Inhibit

AbsCalc
 Ref Trans Schulz reflection method Change Absorption coefieien 133.0 1/cm Thickness 0.2 cm Set 2Theta 38.46 deg. 1/Kt Profile

Defocus file Select Transmission defocus HKL+T
 Defocus(1) functions file C:\CTR\DATA\Aluminum-H-O\AI_random\defocus\DEFOCUS_F.TXT
Make defocus function files by TXT2 Files Normalization degree of a polynomial 0 TenckhoffFitting TextDisp

Defocus(3) function files folder(Calc unbackdefocus) BB185mm Limit Alfa Defocus value Free(LimitValue=0.0)

Defocus(2) function files folder(Calc backdefocus) DSH1.2mm+Schulz+RSH5mm Search minimum EqualAngleRp%(Cubic only) 1/Ra Profile

Smoothing for ADC
 Cycles 2 Weight 10 Disp After connection


Normalization
 CTR Connect Average

CenterData
 Average

OutFiles
 ASC MTexAsc Ras TXT TXT2

ValueODFVF-B ValueODFVF-A **Average** Exit&ODF ODF ALLODFFiles

Filemake success !! CTRHome : C: Select crystal : Cubic 23/08/20

AddingPole.2.19 by CTR PDuser HelperTex CTR

File Help

Files Select(TXT2)
111-7mm_chB00D2S_2.TXT 200-7mm_chB00D2S_2.TXT 220-7mm_chB00D2S_2.TXT 311-7mm_chB00D2S_2.TXT


Holder
C:\A1100P\N0002
C:\A1100P\N0003
C:\A1100P\N0004
C:\A1100P\N0005
C:\A1100P\N0006

Selectnumber: 24 NO2Datapercent 50 % TXT2 Asc

OutFilename
C:\A1100P\N0001\111-7mm_chB00D2S_4Average_2.TXT

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SelectFiles Start **Exit**

ODFPoleFigure2S 4.08 by CTR PDuser HelperTex CTR

File Linear(absolute)3D ToolKit Help InitSet BGMode Measure Condition Free OverlapRevision MinimumMode Rp% Normalization

Files select
ASC(RINT-PC) 111-7mm.ASC 200-7mm.ASC 220-7mm.ASC 311-7mm.ASC

Holder
Holder C:\A1100P\NO004
Rename C:\A1100P\NO005
Renamw-delete C:\A1100P\NO006

Calcration Condition
Previous Next C:\A1100P\NO001\111-7mm.ASC hkl 1,1,1 Change

Background delete mode
 DoubleMo... SingleMode LowMode HighMode Nothing BG defocus DSH1.2mm+Schulz+RSH5mm Minimum mo...
 Minimum(α β) MinimumAverage(α X) 0.5 Trans blinds angle 30.0

Smoothing
 + α 3 Arithmetic mean Disp
RD 0.0 Interporation Full Disp

Peak slit 7.0 mm BG Slit 7.0 mm PeakSlit / BGS... BG Scope 80.0 deg. 90.0 deg. Set SetAll Disp α Inhibit

AbsCalc
 Ref Trans Schulz reflection method Change Absorption coefieien 133.0 1/cm Thickness 0.2 cm Set 2Theta 38.46 deg. 1/Kt Profile

Defocus file Select Transmission defocus HKL+T
 Defocus(1) functions file C:\CTR\DATA\Aluminum-H-O\AI\random\defocus\DEFOCUS_F.TXT
Make defocus function files by TXT2 Files Normalization degree of a polynomial 0 TenckhoffFitting TextDisp

Defocus(3) function files folder(Calc unbackdefocus) BB185mm Limit Alfa Defocus value Free(LimitValue=0.0)

Defocus(2) function files folder(Calc backdefocus) DSH1.2mm+Schulz+RSH5mm Search minimum EqualAngleRp%(Cubic only) 1/Ra Profile

Smoothing for ADC
 Cycles 2 Weight 10 Disp
 After connection

Normalization CenterData
 CTR Connect Average

OutFiles
 ASC MTexAsc Ras TXT TXT2

Cancel Calc Connect

ValueODFVF-B ValueODFVF-A Average **Exit&ODF** ODF ALLODFFiles

Filemake success !! CTRHome : C: Select crystal : Cubic 23/08/20

PfToODF3 8.55 by CTR PDuser HelperTex CTR

File Option Symmetric Software Data Help

Lattice constant
Material Aluminum.txt Initialize Start

Structure Code(Symmetries after Schoenfiles) cif 7 - O (cubic)

a 1.0 <=b 1.0 <=c 1.0 alpha 90.0 beta 90.0 gamm 90.0

PF Holder
C:\A1100P\NO001

PF Data

SelectFile(TXT(b.intens),TXT2(a.b.intens))	h,k,l	2Theta	Alpha scope	AlphaS	AlphaE	Select
111-7mm_chB00D2S_6Average_2.TXT	1,1,1	38.46	0.0->75.0	0.0	75.0	<input checked="" type="checkbox"/>
200-7mm_chB00D2S_6Average_2.TXT	2,0,0	44.7	0.0->75.0	0.0	75.0	<input checked="" type="checkbox"/>
220-7mm_chB00D2S_6Average_2.TXT	2,2,0	65.08	0.0->75.0	0.0	75.0	<input checked="" type="checkbox"/>
311-7mm_chB00D2S_6Average_2.TXT	3,1,1	78.22	0.0->75.0	0.0	75.0	<input checked="" type="checkbox"/>
	2,1,1	0.0		0.0	0.0	<input type="checkbox"/>